

For Research Use Only

MPO Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68797-1-PBS

Basic Information

Catalog Number:

68797-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC130476

GeneID (NCBI):

4353

UNIPROT ID:

P05164

Full Name:

myeloperoxidase

Calculated MW:

745 aa, 84 kDa

Purification Method:

Protein G purification

CloneNo.:

1G1H5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50164-1, MPO Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68797-1-PBS. Detection antibody: 68797-2-PBS. Standard:Eg0511. Range: 0.781-100 ng/mL.