

For Research Use Only

APOBEC3G Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68799-2-PBS

Basic Information

Catalog Number: 68799-2-PBS	GenBank Accession Number: BC024268	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 60489	CloneNo.: 1B1A9
Source: Mouse	UNIPROT ID: Q9HC16	
Isotype: IgG2b	Full Name: apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3G	
Immunogen Catalog Number: AG30583	Calculated MW: 46 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50166-1, APOBEC3G Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68799-1-PBS. Detection antibody: 68799-2-PBS. Standard: Ag30583. Range: 0.391-100 ng/mL.