

For Research Use Only

SMPD3 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68801-2-PBS

Basic Information

Catalog Number:

68801-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG22487

GenBank Accession Number:

BC112238

GeneID (NCBI):

55512

UNIPROT ID:

Q9NY59

Full Name:

sphingomyelin phosphodiesterase 3,
neutral membrane (neutral
sphingomyelinase II)

Calculated MW:

655 aa, 71 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1H3D3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50168-1, SMPD3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68801-1-PBS. Detection antibody: 68801-2-PBS. Standard: Ag22487. Range: 0.781-100 ng/mL.