

For Research Use Only

ENT1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68806-2-PBS

Basic Information

Catalog Number:

68806-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG32016

GenBank Accession Number:

NM_001372327.1

GeneID (NCBI):

2030

UNIPROT ID:

Q99808

Full Name:

solute carrier family 29 (nucleoside transporters), member 1

Calculated MW:

50KD

Purification Method:

Protein G Magarose purification

CloneNo.:

3F5A6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50173-1, ENT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68806-1-PBS. Detection antibody: 68806-2-PBS. Standard: Ag32016. Range: 1.563-100 ng/mL.