

For Research Use Only

IL-1R1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68808-2-PBS



Basic Information

Catalog Number:

68808-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26009

GenBank Accession Number:

BC075063

GeneID (NCBI):

3554

UNIPROT ID:

P14778

Full Name:

interleukin 1 receptor, type I

Calculated MW:

569 aa, 65 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D4B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50175-1, IL-1R1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68808-1-PBS. Detection antibody: 68808-2-PBS. Standard: Ag26009. Range: 0.781-50 ng/mL