

For Research Use Only

Vinculin Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68809-1-PBS

Basic Information

Catalog Number: 68809-1-PBS	GenBank Accession Number: BC039174	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7414	CloneNo.: 4B1D7
Source: Mouse	UNIPROT ID: P18206	
Isotype: IgG2a	Full Name: vinculin	
Immunogen Catalog Number: AG3705	Calculated MW: 1133 aa, 124 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50176-1, Vinculin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68809-1-PBS. Detection antibody: 68809-2-PBS. Standard: ag3705. Range: 0.098-6.25 ng/mL