

For Research Use Only

CHMP4B Monoclonal antibody, PBS Only



Catalog Number: 68811-1-PBS

Basic Information

Catalog Number: 68811-1-PBS	GenBank Accession Number: BC033859	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 128866	CloneNo.: 1D4G10
Source: Mouse	UNIPROT ID: Q9H444	
Isotype: IgG1	Full Name: chromatin modifying protein 4B	
Immunogen Catalog Number: AG4544	Calculated MW: 224 aa, 25 kDa	

Applications

Tested Applications:
WB, IF, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68811-1-Ig (CHMP4B antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68811-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HCT 116 cells using CHMP4B antibody (68811-1-Ig, Clone: 1D4G10) at dilution of 1:2000 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1), CL594-phalloidin (red). This data was developed using the same antibody clone with 68811-1-PBS in a different storage buffer formulation.