

For Research Use Only

# ATF1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68813-2-PBS

## Basic Information

|   |  |  |
|---|--|--|
| <b>Catalog Number:</b><br>68813-2-PBS       | <b>GenBank Accession Number:</b><br>BC029619           | <b>Purification Method:</b><br>Protein A Magarose purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>466                           | <b>CloneNo.:</b><br>2D9C7                                      |
| <b>Source:</b><br>Mouse                     | <b>UNIPROT ID:</b><br>P18846                           |  |
| <b>Isotype:</b><br>IgG3                     | <b>Full Name:</b><br>activating transcription factor 1 |  |
| <b>Immunogen Catalog Number:</b><br>AG34935 | <b>Calculated MW:</b><br>271 aa, 29 kDa                |  |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50187-1, ATF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68813-1-PBS. Detection antibody: 68813-2-PBS. Standard: Ag34935. Range: 0.098-12.5 ng/mL