

For Research Use Only

Centrin 2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68814-2-PBS

Basic Information

Catalog Number: 68814-2-PBS	GenBank Accession Number: BC005334	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1069	CloneNo.: 2B7E4
Source: Mouse	UNIPROT ID: P41208	
Isotype: IgG1	Full Name: centrin, EF-hand protein, 2	
Immunogen Catalog Number: AG8567	Calculated MW: 172 aa, 20 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50188-1, Centrin 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68814-1-PBS. Detection antibody: 68814-2-PBS. Standard:null. Range: null.