

For Research Use Only

VMAT2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68816-1-PBS

Basic Information

Catalog Number: 68816-1-PBS	GenBank Accession Number: BC108928	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6571	CloneNo.: 1E8D1
Source: Mouse	UNIPROT ID: Q05940	
Isotype: IgG1	Full Name: solute carrier family 18 (vesicular monoamine), member 2	
Immunogen Catalog Number: AG14971	Calculated MW: 514 aa, 56 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50190-1, VMAT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68816-1-PBS. Detection antibody: 68816-2-PBS. Standard: Ag14971. Range: 0.391-1.563 ng/mL.