

For Research Use Only

# VMAT2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68816-2-PBS

## Basic Information

Catalog Number:

68816-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14971

GenBank Accession Number:

BC108928

GeneID (NCBI):

6571

UNIPROT ID:

Q05940

Full Name:

solute carrier family 18 (vesicular monoamine), member 2

Calculated MW:

514 aa, 56 kDa

Purification Method:

Protein G purification

CloneNo.:

2E5C12

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50190-1, VMAT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68816-1-PBS. Detection antibody: 68816-2-PBS. Standard: Ag14971. Range: 0.391-1.563 ng/mL.