

For Research Use Only

GBA2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68819-1-PBS

Basic Information

Catalog Number: 68819-1-PBS	GenBank Accession Number: BC011363	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 57704	CloneNo.: 2C1D2
Source: Mouse	UNIPROT ID: Q9HCG7	
Isotype: IgG2b	Full Name: glucosidase, beta (bile acid) 2	
Immunogen Catalog Number: AG10606	Calculated MW: 927 aa, 105 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50194-1, GBA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68819-1-PBS. Detection antibody: 68819-2-PBS. Standard: Ag10606. Range: 0.195-6.25 ng/mL