

For Research Use Only

GBA2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68819-2-PBS

Basic Information

Catalog Number:

68819-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG10606

GenBank Accession Number:

BC011363

GeneID (NCBI):

57704

UNIPROT ID:

Q9HCG7

Full Name:

glucosidase, beta (bile acid) 2

Calculated MW:

927 aa, 105 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D10B1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50194-1, GBA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68819-1-PBS. Detection antibody: 68819-2-PBS. Standard: Ag10606. Range: 0.195-6.25 ng/mL