

For Research Use Only

ODF2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68820-2-PBS

Basic Information

Catalog Number: 68820-2-PBS	GenBank Accession Number: BC010629	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4957	CloneNo.: 3D10H11
Source: Mouse	UNIPROT ID: Q5BJF6	
Isotype: IgG1	Full Name: outer dense fiber of sperm tails 2	
Immunogen Catalog Number: AG33139	Calculated MW: 763 aa, 89 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50196-1, ODF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68820-1-PBS. Detection antibody: 68820-2-PBS. Standard: Ag33139. Range: 3.125-100 ng/mL.