

For Research Use Only

CREBZF Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68821-2-PBS

Basic Information

Catalog Number: 68821-2-PBS	GenBank Accession Number: BC093796	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 58487	CloneNo.: 1F2C10
Source: Mouse	UNIPROT ID: Q9NS37	
Isotype: IgG2b	Full Name: CREB/ATF bZIP transcription factor	
Immunogen Catalog Number: AG13648	Calculated MW: 272 aa, 30 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50198-1, CREBZF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68821-1-PBS. Detection antibody: 68821-2-PBS. Standard:Ag29921. Range: 0.781-100 ng/mL.