

For Research Use Only

REG3A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68822-1-PBS

Basic Information

Catalog Number: 68822-1-PBS	GenBank Accession Number: BC036776	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5068	CloneNo.: 1A1A2
Source: Mouse	UNIPROT ID: Q06141	
Isotype: IgG2a	Full Name: regenerating islet-derived 3 alpha	
Immunogen Catalog Number: AG4725	Calculated MW: 175 aa, 19 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50201-1, REG3A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68822-1-PBS. Detection antibody: 68822-2-PBS. Standard: Ag4725. Range: 0.098-100 ng/mL.