

For Research Use Only

DDAH1 Monoclonal antibody, PBS Only



Catalog Number: 68823-1-PBS

Basic Information

Catalog Number:

68823-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG32840

GenBank Accession Number:

NM_012137

GeneID (NCBI):

23576

UNIPROT ID:

O94760

Full Name:

dimethylarginine

dimethylaminohydrolase 1

Calculated MW:

31kd

Purification Method:

Protein G purification

CloneNo.:

1A1A6

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68823-1-Ig (DDAH1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68823-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 68823-1-Ig (DDAH1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68823-1-PBS in a different storage buffer formulation.