

For Research Use Only

# Histone H3.1 Monoclonal antibody, PBS Only



Catalog Number: 68825-1-PBS

## Basic Information

Catalog Number: 68825-1-PBS	GenBank Accession Number: BC066245	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8350	CloneNo.: 2E6D11
Source: Mouse	UNIPROT ID: P68431	
Isotype: IgG2a	Full Name: histone cluster 1, H3a	

## Applications

Tested Applications:  
WB, ELISA  
Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68825-1-Ig (Histone H3.1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68825-1-PBS in a different storage buffer formulation.