

For Research Use Only

CD69 Monoclonal antibody, PBS Only



Catalog Number: 68827-1-PBS

Basic Information

Catalog Number:

68827-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG30963

GenBank Accession Number:

BC007037

GeneID (NCBI):

969

UNIPROT ID:

Q07108

Full Name:

CD69 molecule

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

3D10A2

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68827-1-Ig (CD69 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68827-1-PBS in a different storage buffer formulation.