

For Research Use Only

# Fas/CD95 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68828-3-PBS

## Basic Information

Catalog Number:

68828-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_000043.6

GeneID (NCBI):

355

UNIPROT ID:

P25445-1

Full Name:

Fas (TNF receptor superfamily, member 6)

Calculated MW:

38 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E3D7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50205-2, Fas/CD95 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68828-1-PBS. Detection antibody: 68828-3-PBS. Standard:Eg0877. Range: 1.563-100 ng/mL