

For Research Use Only

KLHL21 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68830-1-PBS

Basic Information

Catalog Number:

68830-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10845

GenBank Accession Number:

BC034039

GeneID (NCBI):

9903

UNIPROT ID:

Q9UJP4

Full Name:

kelch-like 21 (Drosophila)

Calculated MW:

597 aa, 67 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E4H2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50207-1, KLHL21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68830-1-PBS. Detection antibody: 68830-2-PBS. Standard:Ag10845. Range: 0.781-100 ng/mL.