

For Research Use Only

KLHL21 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68830-2-PBS

Basic Information

Catalog Number: 68830-2-PBS	GenBank Accession Number: BC034039	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9903	CloneNo.: 1F12G5
Source: Mouse	UNIPROT ID: Q9UJP4	
Isotype: IgG1	Full Name: kelch-like 21 (Drosophila)	
Immunogen Catalog Number: AG10845	Calculated MW: 597 aa, 67 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50207-1, KLHL21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68830-1-PBS. Detection antibody: 68830-2-PBS. Standard:Ag10845. Range: 0.781-100 ng/mL.