

For Research Use Only

# C16orf63 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68832-1-PBS

## Basic Information

Catalog Number:

68832-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24071

GenBank Accession Number:

BC022321

GeneID (NCBI):

123811

UNIPROT ID:

Q96NB1

Full Name:

chromosome 16 open reading frame

63

Purification Method:

Protein G Magarose purification

CloneNo.:

1F1B7

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50210-1, C16orf63 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68832-1-PBS. Detection antibody: HRP-conjugated 68832-2-PBS. Standard: Ag24071. Range: 0.049-50 ng/mL.