

For Research Use Only

CALM1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68833-1-PBS

Basic Information

Catalog Number:

68833-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG35001

GenBank Accession Number:

BC000454

GeneID (NCBI):

801

UNIPROT ID:

PODP23

Full Name:

calmodulin 1 (phosphorylase kinase, delta)

Calculated MW:

17 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2H9B7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50211-1, CALM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68833-1-PBS. Detection antibody: HRP-conjugated 68833-2-PBS. Standard: Ag35001. Range: 0.049-50 ng/mL.