

For Research Use Only

IGFBP4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68834-3-PBS

Basic Information

Catalog Number: 68834-3-PBS	GenBank Accession Number: BC016041	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3487	CloneNo.: 1G3C5
Source: Mouse	UNIPROT ID: P22692	
Isotype: IgG1	Full Name: insulin-like growth factor binding protein 4	
Immunogen Catalog Number: AG34800	Calculated MW: 28 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50212-2, IGFBP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68834-3-PBS. Detection antibody: HRP-conjugated 68834-4-PBS. Standard: Ag34800. Range: 0.049-50 ng/mL.