

For Research Use Only

# SHBG Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68835-2-PBS

## Basic Information

Catalog Number:

68835-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC101785

GeneID (NCBI):

6462

UNIPROT ID:

P04278

Full Name:

sex hormone-binding globulin

Calculated MW:

402 aa, 44 kDa

Purification Method:

Protein G purification

CloneNo.:

1B9E11

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50219-1, SHBG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68835-2-PBS. Detection antibody: 68835-3-PBS. Standard: Eg0369. Range: 0.625-80 ng/mL