

For Research Use Only

# NF2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68839-1-PBS

## Basic Information

<b>Catalog Number:</b> 68839-1-PBS	<b>GenBank Accession Number:</b> BC020257	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 4771	<b>CloneNo.:</b> 3D2A2
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P35240	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> neurofibromin 2 (merlin)	
<b>Immunogen Catalog Number:</b> AG16265	<b>Calculated MW:</b> 595 aa, 70 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50217-1, NF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68839-1-PBS. Detection antibody: 68839-2-PBS. Standard: Ag16265. Range: 0.098-100 ng/mL.