

For Research Use Only

NF2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68839-5-PBS

Basic Information

Catalog Number: 68839-5-PBS	GenBank Accession Number: BC020257	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4771	CloneNo.: 3A11E5
Source: Mouse	UNIPROT ID: P35240	
Isotype: IgG1	Full Name: neurofibromin 2 (merlin)	
Immunogen Catalog Number: AG16265	Calculated MW: 595 aa, 70 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50217-3, NF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68839-5-PBS. Detection antibody: 68839-6-PBS. Standard: Ag16265. Range: 0.391-100 ng/mL.