

For Research Use Only

# VPREB3 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68841-2-PBS

## Basic Information

<b>Catalog Number:</b> 68841-2-PBS	<b>GenBank Accession Number:</b> BC020666	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 29802	<b>CloneNo.:</b> 2H9C6
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9UKI3	
<b>Isotype:</b> IgG1	<b>Full Name:</b> pre-B Lymphocyte 3	
<b>Immunogen Catalog Number:</b> AG2528	<b>Calculated MW:</b> 123 aa, 14 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50225-1, VPRESB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68841-1-PBS. Detection antibody: 68841-2-PBS. Standard: Ag2528. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50225-2, VPRESB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68841-3-PBS. Detection antibody: 68841-2-PBS. Standard: Ag2528. Range: 0.098-100 ng/mL.