

For Research Use Only

PROS1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68843-2-PBS

Basic Information

Catalog Number:

68843-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10503

GenBank Accession Number:

BC015801

GeneID (NCBI):

5627

UNIPROT ID:

P07225

Full Name:

protein S (alpha)

Calculated MW:

75 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D4E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50227-1, PROS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68843-1-PBS. Detection antibody: 68843-2-PBS. Standard: Ag10503. Range: 1.563-100 ng/mL.