

For Research Use Only

PINK1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68848-2-PBS

Basic Information

Catalog Number:

68848-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG19825

GenBank Accession Number:

BC028215

GeneID (NCBI):

65018

UNIPROT ID:

Q9BXM7

Full Name:

PTEN induced putative kinase 1

Calculated MW:

581 aa, 63 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2F9G12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50236-1, PINK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68848-1-PBS. Detection antibody: 68848-2-PBS. Standard: Ag19825. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50236-2, PINK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68848-3-PBS. Detection antibody: 68848-2-PBS. Standard: Ag19825. Range: 0.391-100 ng/mL.