

For Research Use Only

PINK1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68848-3-PBS



Basic Information

Catalog Number: 68848-3-PBS	GenBank Accession Number: BC028215	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 65018	CloneNo.: 4D1C1
Source: Mouse	UNIPROT ID: Q9BXM7	
Isotype: IgG2a	Full Name: PTEN induced putative kinase 1	
Immunogen Catalog Number: AG19825	Calculated MW: 581 aa, 63 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50236-2, PINK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68848-3-PBS. Detection antibody: 68848-2-PBS. Standard: Ag19825. Range: 0.391-100 ng/mL.