

For Research Use Only

ADAR1 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 68850-2-PBS

Basic Information

Catalog Number: 68850-2-PBS	GenBank Accession Number: BC038227	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 103	CloneNo.: 3A1C7
Source: Mouse	UNIPROT ID: P55265	
Isotype: IgG1	Full Name: adenosine deaminase, RNA-specific	
Immunogen Catalog Number: AG5609	Calculated MW: 136 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50240-1, ADAR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68850-1-PBS. Detection antibody: 68850-2-PBS. Standard:Ag5609. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50240-2, ADAR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68850-2-PBS. Detection antibody: 68850-3-PBS. Standard:Ag5609. Range: 0.781-100 ng/mL