

For Research Use Only

CHGA Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68851-1-PBS

Basic Information

Catalog Number: 68851-1-PBS	GenBank Accession Number: BC006459	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1113	CloneNo.: 1D11G12
Source: Mouse	UNIPROT ID: P10645	
Isotype: IgG2a	Full Name: chromogranin A (parathyroid secretory protein 1)	
Immunogen Catalog Number: AG19780	Calculated MW: 51 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50241-1, CHGA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68851-1-PBS. Detection antibody: 60135-3-PBS. Standard:Ag19780. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50241-3, CHGA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68851-1-PBS. Detection antibody: 60135-2-PBS. Standard:Ag19780. Range: 1.563-100 ng/mL.