

For Research Use Only

TSPAN3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68852-3-PBS

Basic Information

Catalog Number: 68852-3-PBS	GenBank Accession Number: BC011206	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10099	CloneNo.: 1A10E3
Source: Mouse	UNIPROT ID: O60637	
Isotype: IgG1	Full Name: tetraspanin 3	
Immunogen Catalog Number: AG3689	Calculated MW: 253 aa, 28 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50243-2, TSPAN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68852-3-PBS. Detection antibody: 68852-2-PBS. Standard: Ag3689. Range: 0.391-6.25 ng/mL



Cytometric bead array standard curve of MP50243-3, TSPAN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68852-3-PBS. Detection antibody: 68852-4-PBS. Standard: Ag3689. Range: 0.098-3.125 ng/mL