

For Research Use Only

CLCC1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68853-1-PBS

Basic Information

| | | |
|---|---|--|
| Catalog Number: 68853-1-PBS | GenBank Accession Number: BC002939 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 23155 | CloneNo.: 3F11E12 |
| Source: Mouse | UNIPROT ID: Q96566 | |
| Isotype: IgG1 | Full Name: chloride channel CLIC-like 1 | |
| Immunogen Catalog Number: AG24682 | Calculated MW: 551 aa, 62 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50244-1, CLCC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68853-1-PBS. Detection antibody: 68853-2-PBS. Standard: Ag24682. Range: 1.563-100 ng/mL.