

For Research Use Only

# CLCC1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68853-2-PBS

## Basic Information

<b>Catalog Number:</b> 68853-2-PBS	<b>GenBank Accession Number:</b> BC002939	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 23155	<b>CloneNo.:</b> 2G5A11
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q96566	
<b>Isotype:</b> IgG3	<b>Full Name:</b> chloride channel CLIC-like 1	
<b>Immunogen Catalog Number:</b> AG24682	<b>Calculated MW:</b> 551 aa, 62 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50244-1, CLCC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68853-1-PBS. Detection antibody: 68853-2-PBS. Standard: Ag24682. Range: 1.563-100 ng/mL.