

For Research Use Only

CLCC1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68853-6-PBS

Basic Information

Catalog Number:

68853-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG24682

GenBank Accession Number:

BC002939

GeneID (NCBI):

23155

UNIPROT ID:

Q96566

Full Name:

chloride channel CLIC-like 1

Calculated MW:

551 aa, 62 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3F10C9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50244-4, CLCC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68853-3-PBS. Detection antibody: 68853-6-PBS. Standard:Ag24682. Range: 0.781-100 ng/mL.