

For Research Use Only

RPL22 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68854-2-PBS

Basic Information

| | | |
|--------------------------------------|---------------------------------------|--|
| Catalog Number: 68854-2-PBS | GenBank Accession Number: BC058887 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 6146 | CloneNo.: 3B7H4 |
| Source: Mouse | UNIPROT ID: P35268 | |
| Isotype: IgG3 | Full Name: ribosomal protein L22 | |
| Immunogen Catalog Number: AG30115 | Calculated MW: 15 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50247-1, RPL22 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68854-1-PBS. Detection antibody: 68854-2-PBS. Standard: Ag30115. Range: 1.563-50 ng/mL