

For Research Use Only

# CLSTN3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68855-1-PBS

## Basic Information

Catalog Number:

68855-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14897

GenBank Accession Number:

BC039075

GeneID (NCBI):

9746

UNIPROT ID:

Q9BQT9

Full Name:

calcsyntenin 3

Calculated MW:

968 aa, 107 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D8E4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50250-1, CLSTN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68855-1-PBS. Detection antibody: 68855-2-PBS. Standard: Ag14897. Range: 0.391-25 ng/mL