

For Research Use Only

ER Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68857-3-PBS

Basic Information

Catalog Number: 68857-3-PBS	GenBank Accession Number: BC128573	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2099	CloneNo.: 1D7A11
Source: Mouse	UNIPROT ID: P03372	
Isotype: IgG1	Full Name: estrogen receptor 1	
Immunogen Catalog Number: AG28104	Calculated MW: 595 aa, 66 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50252-2, ER Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68857-3-PBS. Detection antibody: 68857-4-PBS. Standard:Ag28104. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50252-3, ER Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68857-3-PBS. Detection antibody: 68857-5-PBS. Standard:Ag28104. Range: 3.125-100 ng/mL.