

For Research Use Only

ER Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68858-4-PBS

Basic Information

Catalog Number: 68858-4-PBS	GenBank Accession Number: BC128573	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2099	CloneNo.: 4H6D6
Source: Mouse	UNIPROT ID: P03372	
Isotype: IgG1	Full Name: estrogen receptor 1	
Immunogen Catalog Number: AG15738	Calculated MW: 595 aa, 66 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50252-5, ER Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68858-3-PBS. Detection antibody: 68858-4-PBS. Standard:Ag15738. Range: 1.563-100 ng/mL.