

For Research Use Only

ER Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68858-6-PBS

Basic Information

Catalog Number:

68858-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15738

GenBank Accession Number:

BC128573

GeneID (NCBI):

2099

UNIPROT ID:

P03372

Full Name:

estrogen receptor 1

Calculated MW:

595 aa, 66 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G3B3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50252-6, ER Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68858-5-PBS. Detection antibody: 68858-6-PBS. Standard: Ag15738. Range: 1.563-100 ng/mL.