

For Research Use Only

# MRP3/ABCC3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68861-1-PBS



## Basic Information

Catalog Number:

68861-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19801

GenBank Accession Number:

BC137348

GeneID (NCBI):

8714

UNIPROT ID:

O15438

Full Name:

ATP-binding cassette, sub-family C (CFTR/MRP), member 3

Calculated MW:

1527 aa, 169 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B8G6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50255-1, ABCC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68861-1-PBS. Detection antibody: 68861-2-PBS. Standard: Ag19801. Range: 0.781-100 ng/mL.