

For Research Use Only

CHMP1A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68862-1-PBS

Basic Information

Catalog Number: 68862-1-PBS	GenBank Accession Number: BC007527	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5119	CloneNo.: 1H9D10
Source: Mouse	UNIPROT ID: Q9HD42	
Isotype: IgG3	Full Name: chromatin modifying protein 1A	
Immunogen Catalog Number: AG8376	Calculated MW: 196 aa, 22 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50256-1, CHMP1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68862-1-PBS. Detection antibody: 68862-2-PBS. Standard: Ag8376. Range: 0.781-100 ng/mL.