

For Research Use Only

# CHMP1A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68862-4-PBS

## Basic Information

|  |   |  |
|--|---|--|
| <b>Catalog Number:</b><br>68862-4-PBS      | <b>GenBank Accession Number:</b><br>BC007527        | <b>Purification Method:</b><br>Protein G Magarose purification |
| <b>Size:</b><br>1 mg/ml                    | <b>GeneID (NCBI):</b><br>5119                       | <b>CloneNo.:</b><br>1C8D4                                      |
| <b>Source:</b><br>Mouse                    | <b>UNIPROT ID:</b><br>Q9HD42                        |  |
| <b>Isotype:</b><br>IgG1                    | <b>Full Name:</b><br>chromatin modifying protein 1A |  |
| <b>Immunogen Catalog Number:</b><br>AG8376 | <b>Calculated MW:</b><br>196 aa, 22 kDa             |  |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50256-2, CHMP1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68862-3-PBS. Detection antibody: 68862-4-PBS. Standard:Ag8376. Range: 0.781-100 ng/mL.