

For Research Use Only

PLK4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68864-1-PBS

Basic Information

Catalog Number: 68864-1-PBS	GenBank Accession Number: BC036023	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10733	CloneNo.: 2G9B7
Source: Mouse	UNIPROT ID: O00444	
Isotype: IgG1	Full Name: polo-like kinase 4 (Drosophila)	
Immunogen Catalog Number: AG21287	Calculated MW: 970 aa, 109 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50258-1, PLK4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68864-1-PBS. Detection antibody: 68864-2-PBS. Standard: Ag21287. Range: 0.391-50 ng/mL