

For Research Use Only

LSDP5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68865-2-PBS

Basic Information

Catalog Number: 68865-2-PBS	GenBank Accession Number: BC131524	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 440503	CloneNo.: 1E5G12
Source: Mouse	UNIPROT ID: Q00G26	
Isotype: IgG1	Full Name: lipid storage droplet protein 5	
Immunogen Catalog Number: AG35548	Calculated MW: 463 aa, 51 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50260-1, LSDPS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68865-1-PBS. Detection antibody: 68865-2-PBS. Standard: Ag35548. Range: 1.563-100 ng/mL.