

For Research Use Only

TIMM8A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68866-2-PBS

Basic Information

Catalog Number:

68866-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33190

GenBank Accession Number:

BC015093

GeneID (NCBI):

1678

UNIPROT ID:

O60220

Full Name:

translocase of inner mitochondrial membrane 8 homolog A (yeast)

Calculated MW:

11 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H3E2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50261-1, TIMM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68866-1-PBS. Detection antibody: 68866-2-PBS. Standard: Ag33190. Range: 1.563-100 ng/mL.