

For Research Use Only

# TIMM8A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68866-4-PBS

## Basic Information

Catalog Number:

68866-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33190

GenBank Accession Number:

BC015093

GeneID (NCBI):

1678

UNIPROT ID:

O60220

Full Name:

translocase of inner mitochondrial membrane 8 homolog A (yeast)

Calculated MW:

11 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B8E5

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50261-3, TIMM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68866-1-PBS. Detection antibody: 68866-4-PBS. Standard: Ag33190. Range: 1.563-100 ng/mL.